

**Bottiglia con tappo a vite, 200 ml, in PC-CB****Herolab**

Bottiglia con tappo a vite, 200 ml, in PC-CB. Dim. diam. 56 x 112 mm. Per centrifuga. Rpm: 10.000.



Code	Description	Packaging
LLG06522461	Bottiglia con tappo a vite, 200 ml, in PC-CB	6 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

