



Inserto di ricambio in PTFE, per tappo a vite, diam. 13 mm, conf. 20

DWK Life Sciences

Inserto di ricambio in PTFE, per tappo a vite, diam. 13 mm.



Code	Description	Packaging
LLG06524332	Pyrex Screwcap Spare Inserts PTFE, 13 mm pack of 20	20 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

