



Igrometro Wg. blu, con termometro 0...+35 °C, lungo circa 330 mm, range 1,000-1,200/ div. 0,002 g/ml, temp. di riferimento 20°C

GECO Gering

Igrometro Wg. blu, con termometro 0...+35 °C, lungo circa 330 mm, range 1,000-1,200/ div. 0,002 g/ml, temp. di riferimento 20°C.



Code	Description	Packaging
LLG07970837	Igrometro con termometro Wg. blu, range 0...+35 °C, lungo circa 330 mm, 1,000-1,200 in 0,002 g/ml, Tp. 20°C	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

