



SCB N9-H, bl, Sil bg/PTFE w(-), 45°, 1.3

Macherey Nagel

N 9 PP screw cap (bonded), blue, center hole Silicone beige/PTFE white, slit Hardness: 45° shore A Thickness: 1.3 mm pack of 100



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL702027 | SCB N9-H, bl, Sil bg/PTFE w(-), 45°, 1.3 | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

