



SC N13-nH, blk, Sil w/PTFE r, 40°, 1.3

Macherey Nagel

N 13 PP screw cap, black, closed top Silicone white/PTFE red Hardness: 40°  
shore A Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702052	SC N13-nH, blk, Sil w/PTFE r, 40°, 1.3	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

