



SCB N24-nH, w, Sil w/PTFE bg, 45°, 3.2

Macherey Nagel

N 24 PP screw cap (bonded), white, closed top Silicone white/PTFE beige
Hardness: 45° shore A Thickness: 3.2 mm pack of 100



Code	Description	Packaging
MNL702059	SCB N24-nH, w, Sil w/PTFE bg, 45°, 3.2	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

