



SCB N24-H, w, RR or/TEF, 45°, 2.5

Macherey Nagel

N 24 PP screw cap (bonded), white, center hole Red Rubber/TEF colourless Hardness: 45°
shore A Thickness: 2.5 mm pack of 100



Code	Description	Packaging
MNL702073	SCB N24-H, w, RR or/TEF, 45°, 2.5	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

