



CC N20-H, si, Sil w/PTFE bg, 40°, 3.0

Macherey Nagel

N 20 Aluminium crimp cap, silver, center hole Silicone white / PTFE beige
Hardness: 40° shore A Thickness: 3 mm pack of 100



Code	Description	Packaging
MNL702094	CC N20-H, si, Sil w/PTFE bg, 40°, 3.0	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

