



**SCB N18-nH, w, Sil w/PTFE bl, 55°, 1.3**

**Macherey Nagel**

N 18 PP screw cap (bonded, 18-400), white, closed top Silicone white/PTFE  
blue Hardness: 55° shore A Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702115	SCB N18-nH, w, Sil w/PTFE bl, 55°, 1.3	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

