



SCB N20-nH, w, Sil w/PTFE bl, 55°,1.3

Macherey Nagel

N 20 PP screw cap (bonded, 20-400), white, closed top Silicone white/PTFE
blue Hardness: 55° shore A Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702116	SCB N20-nH, w, Sil w/PTFE bl, 55°,1.3	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

