



PE-Cap N20-DIN-H, tr, NR or/TEF,45°,1.3

Macherey Nagel

N 20 PE Cap for crimp neck vials with flat DIN crimp neck N 20 transparent,
center hole Natural Rubber red-orange/TEF colourless Hardness: 45° shore A
Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702129	PE-Cap N20-DIN-H, tr, NR or/TEF,45°,1.3	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

