



Mg SC N18-H,si, Sil w/PTFE bl(-) 45°,1.5

Macherey Nagel

N 18 magnetic screw cap, silver, center hole Silicone white/PTFE blue, slit Hardness: 45° shore
A Thickness: 1.5 mm pack of 100



Code	Description	Packaging
MNL702136	Mg SC N18-H,si, Sil w/PTFE bl(-) 45°,1.5	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

