



**Mg SC N18-nH,si, Sil w/PTFE bl, 45°, 1.5**

**Macherey Nagel**

N 18 magnetic screw cap, silver, closed top Silicone white/PTFE blue Hardness:  
45° shore A Thickness: 1.5 mm pack of 100



Code	Description	Packaging
MNL702140	Mg SC N18-nH,si, Sil w/PTFE bl, 45°, 1.5	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

