



Mg SC N18-nH,si, Sil w/PTFE bl, 45°, 1.5

Macherey Nagel

N 18 magnetic screw cap, silver, closed top Silicone white/PTFE blue Hardness: 45° shore A
Thickness: 1.5 mm pack of 100



Code	Description	Packaging
MNL702140	Mg SC N18-nH,si, Sil w/PTFE bl, 45°, 1.5	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

