



CC N20-H, si, Sil w/FEP/Alu si, 45°, 3.2

Macherey Nagel

N 20 Aluminium crimp cap, silver, center hole Silicone white / FEP-/Aluminium foil silver
Hardness: 45° shore A Thickness: 3.2 mm pack of 100



Code	Description	Packaging
MNL702145	CC N20-H, si, Sil w/FEP/Alu si, 45°, 3.2	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

