



SC N9-H, mag, si, Sil w/PTFE bl,(-)55°,1

Macherey Nagel

N 9 magnetic screw cap, silver, center hole Silicone white/PTFE blue, slit Hardness: 55° shore A
Thickness: 1.0 mm pack of 100



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL702156 | SC N9-H, mag, si, Sil w/PTFE bl,(-)55°,1 | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

