



**SCB N24-H/protcap,w,Sil nat/PTFE,45°,3.2**

**Macherey Nagel**

N 24 PP screw cap (bonded), white, 15 mm center hole, with removable protection lid Silicone natural/PTFE colourless Hardness: 45° shore A Thickness: 3.2 mm pack of 100



Code	Description	Packaging
MNL702168	SCB N24-H/protcap,w,Sil nat/PTFE,45°,3.2	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

