



CC N20-H,si,But dgr/ct.PTFEgr,50°,3

Macherey Nagel

N 20 Aluminium crimp cap, silver, center hole Butyl dark gray/central PTFE-lamination gray
Hardness: 50° shore A Thickness: 3 mm pack of 100



Code	Description	Packaging
MNL70234	CC N20-H,si,But dgr/ct.PTFEgr,50°,3	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

