



CC N20-H,si, Sil bltr/ct. PTFE, 45°, 3.0

Macherey Nagel

N 20 Aluminium crimp cap, silver, center hole Silicone blue transparent/central PTFE-lamination colorless Hardness: 45° shore A Thickness: 3 mm pack of 100



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL702340 | CC N20-H,si, Sil bltr/ct. PTFE, 45°, 3.0 | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

