



PE-Cap N20-DIN-H tr, Bu/PTFE 55° 1.3 na

Macherey Nagel

N 20 PE Cap for crimp neck vials with flat DIN crimp neck N 20 transparent, center hole Butyl beige/PTFE gray Hardness: 55° shore A Thickness: 1.3 mm (unassembled) pack of 100



| Code | Description | Packaging |
|----------|---|-----------|
| MNL70240 | PE-Cap N20-DIN-H tr, Bu/PTFE 55° 1.3 na | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

