



SR N11-H, Ibl, s, Sil w/Polyimide or,1.0

Macherey Nagel

N 11 PE snap ring cap, soft, light blue, center hole Silicone white/Polyimide orange Hardness: 50° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL702403	SR N11-H, lbl, s, Sil w/Polyimide or,1.0	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

