



Tappo in PE N20-bevT-H tr, Bu/PTFE 55° 1.3 na

Macherey Nagel

N 20 PE Cap for crimp neck vials with beveled top N 20 transparent, center hole  
Butyl beige/PTFE gray Hardness: 55° shore A Thickness: 1.3 mm  
(unassembled) pack of 100



Code	Description	Packaging
MNL70242	PE-Cap N20-bevT-H tr, Bu/PTFE 55° 1.3 na	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

