



**Setto alte temperature N 17 HT, rosso**

**Macherey Nagel**

High temperature septa N 17 HT OD: 17 mm, thickness: 3 mm, pack of 50



| Code      | Description                         | Packaging |
|-----------|-------------------------------------|-----------|
| MNL702632 | High temperature septa N 17 HT, red | 50 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

