



CenterGuide BTO Septa Shimadzu

Macherey Nagel

CenterGuide BTO Septa Shimadzu GC septa, BTO OD: Shimadzu plug, thickness: Shimadzu plug max. temp. 300°C, coated blister pack of 25



| Code | Description | Packaging |
|-----------|--------------------------------|-----------|
| MNL702649 | CenterGuide BTO Septa Shimadzu | 25 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

