



CenterGuide AG3 Septa N 11.5

Macherey Nagel

CenterGuide AG3 Septa N 11.5 GC septa, Advanced Green 3 OD: 11.5 mm, thickness: 3 mm
max. temp. 400°C, coated blister pack of 25



| Code | Description | Packaging |
|-----------|------------------------------|-----------|
| MNL702658 | CenterGuide AG3 Septa N 11.5 | 25 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

