



**CenterGuide AG3 Septa Shimadzu**

**Macherey Nagel**

CenterGuide AG3 Septa Shimadzu OD: GC Septa, Advanced Green 3 OD:  
Shimadzu plug, thickness: Shimadzu plug max. temp. 300°C, coated blister pack  
of 25



Code	Description	Packaging
MNL702659	CenterGuide AG3 Septa Shimadzu	25 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

