



CenterGuide AG3 Septa Shimadzu

Macherey Nagel

CenterGuide AG3 Septa Shimadzu OD: GC Septa, Advanced Green 3 OD: Shimadzu plug,
thickness: Shimadzu plug max. temp. 300°C, coated blister pack of 25



Code	Description	Packaging
MNL702659	CenterGuide AG3 Septa Shimadzu	25 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

