



Setti per porta di iniezione GC, qualità Premium HT, Marathon, CenterGuide, 9 mm

Macherey Nagel

CenterGuide Marathon Septa N 9 GC septa, Marathon OD: 9 mm, thickness: 3 mm max. temp. 400°C, coated blister pack of 25



Code	Description	Packaging
MNL702660	Setti per porta di iniezione GC, qualità Premium HT, Marathon, CenterGuide, 9 mm	25 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

