



CenterGuide Marathon Septa N 11

Macherey Nagel

CenterGuide Marathon Septa N 11 GC septa, Marathon OD: 11 mm, thickness:
3 mm max. temp. 400°C, coated blister pack of 25



Code	Description	Packaging
MNL702661	CenterGuide Marathon Septa N 11	25 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

