



**CenterGuide Marathon Septa N 11**

**Macherey Nagel**

CenterGuide Marathon Septa N 11 GC septa, Marathon OD: 11 mm, thickness: 3 mm max.  
temp. 400°C, coated blister pack of 25



Code	Description	Packaging
MNL702661	CenterGuide Marathon Septa N 11	25 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

