



SR N11-H, lbl, s, Sil w/PTFE r, 45°, 1.0

Macherey Nagel

N 11 PE snap ring cap, soft, light blue, center hole Silicone white/PTFE red
Hardness: 45° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL7027102080	SR N11-H, lbl, s, Sil w/PTFE r, 45°, 1.0	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

