



SR N11-H, bl, Sil w/PTFE bl (+), 55° 1.0

Macherey Nagel

N 11 PE snap ring cap, blue, center hole Silicone white/PTFE blue, cross-slit
Hardness: 55° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL7027172	SR N11-H, bl, Sil w/PTFE bl (+), 55° 1.0	100 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

