



**SR N11-H, lbl,s,Sil w/PTFE bl (+),55° 1.**

**Macherey Nagel**

N 11 PE snap ring cap, soft, light blue, center hole Silicone white/PTFE blue, cross-slit Hardness: 55° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL7027172080	SR N11-H, lbl,s,Sil w/PTFE bl (+),55° 1.	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

