



Setto N20 Sil w/PTFE bg, 45°, 3.0

Macherey Nagel

N 20 septa Silicone white/PTFE beige Hardness: 45° shore A Thickness: 3 mm
pack of 100



Code	Description	Packaging
MNL70278	Septa N20 Sil w/PTFE bg, 45°, 3.0	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

