



CC N8-H si, PTFEr/Sil w/PTFEr, 50°, 1.0

Macherey Nagel

N 8 Aluminium crimp cap, silver, center hole PTFE red/Silicone white/PTFE red
Hardness: 50° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL702878	CC N8-H si, PTFEr/Sil w/PTFEr, 50°, 1.0	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

