



Mag CC N11-H, go, Sil w/PTFE r, 45°, 1.3

Macherey Nagel

N 11 Magnetic crimp cap, gold, center hole Silicone white/PTFE red Hardness:
45° shore A Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702879	Mag CC N11-H, go, Sil w/PTFE r, 45°, 1.3	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

