



SC N13-H, blk, Sil w/PTFE r, 40°, 1.3

Macherey Nagel

N 13 PP screw cap, black, center hole Silicone white/PTFE red Hardness: 40° shore A
Thickness: 1.3 mm pack of 100



Code	Description	Packaging
MNL702926	SC N13-H, blk, Sil w/PTFE r, 40°, 1.3	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

