



MCC N20-H, si, But lgr/PTFE dgr, 50°,3.0

Macherey Nagel

N 20 Magnetic crimp cap, silver, 8 mm center hole Butyl light gray/PTFE dark gray Hardness: 50° shore A Thickness: 3 mm pack of 100



Code	Description	Packaging
MNL702928	MCC N20-H, si, But lgr/PTFE dgr, 50°,3.0	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

