



**MCC N20-H, si, Sil bltr/PTFE, 45°, 3.0**

**Macherey Nagel**

N 20 Magnetic crimp cap, silver, 8 mm center hole Silicone blue transparent/PTFE colourless  
Hardness: 45° shore A Thickness: 3 mm pack of 100



Code	Description	Packaging
MNL702929	MCC N20-H, si, Sil bltr/PTFE, 45°, 3.0	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

