



Inserti 0.3, c, 6x40, metal spring req.

Macherey Nagel

0.3 mL Insert outer diameter: 6 mm, outer height: 40 mm clear, conical, metal spring 702974
necessary pack of 100



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL702972 | Insert 0.3, c, 6x40, metal spring req. | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

