



CC N11-H, si, PTFEr/Sil w/PTFEr, 50° 1.0

Macherey Nagel

N 11 Aluminium crimp cap, silver, center hole PTFE red/Silicone white/PTFE red
Hardness: 50° shore A Thickness: 1.0 mm pack of 100



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL702995 | CC N11-H, si, PTFEr/Sil w/PTFEr, 50° 1.0 | 100 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

