



**CC N11-H, bl, PTFEr/Sil w/PTFEr, 50° 1.0**

**Macherey Nagel**

N 11 Aluminium crimp cap, blue, center hole PTFE red/Silicone white/PTFE red  
Hardness: 50° shore A Thickness: 1.0 mm pack of 100



Code	Description	Packaging
MNL7029953	CC N11-H, bl, PTFEr/Sil w/PTFEr, 50° 1.0	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

