



Flash-Chromat.column 450x40 mm

Macherey Nagel

Flash chromatography column, glass complete with adaptor, PTFE valve length:
450 mm, ID: 40 mm



Code	Description	Packaging
MNL727407	Flash-Chromat.column 450x40 mm	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

