



**Flash-Chromat.column 450x40 mm**

**Macherey Nagel**

Flash chromatography column, glass complete with adaptor, PTFE valve length: 450 mm, ID: 40 mm



Code	Description	Packaging
MNL727407	Flash-Chromat.column 450x40 mm	1 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

