

**Flash-Chromatography kit II, 450x40****Macherey Nagel**

Flash chromatography kit II consists of: glass column 450/40 2 l eluent reservoir, 100 g silica gel 60, 1 m PTFE tubing, 500 g seasand, glass fiber wadding, silanized



Code	Description	Packaging
MNL727451	Flash-Chromatography kit II, 450x40	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

