



**CHROMAFIL Xtra GF/PA-20/25**

**Macherey Nagel**

CHROMAFIL Xtra disposable syringe filters GF/PA-20/25 membrane material: glass fiber/polyamide pore dia.: 1.0/0.2  $\mu\text{m}$ , membrane dia.: 25 mm colorless PP case, labeled pack of 100



Code	Description	Packaging
MNL729260	CHROMAFIL Xtra GF/PA-20/25	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

