

**CHROMAFIL Xtra GF/PTFE-45/25****Macherey Nagel**

CHROMAFIL Xtra disposable syringe filters GF/PTFE-45/25 membrane material:  
glass fiber/polytetrafluoroethylene pore dia.: 1.0/0.45  $\mu\text{m}$ , membrane dia.: 25 mm  
colorless PP case, labeled pack of 100



Code	Description	Packaging
MNL729271	CHROMAFIL Xtra GF/PTFE-45/25	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

