

**CHROMAFIL Xtra GF/PTFE-45/25****Macherey Nagel**

CHROMAFIL Xtra disposable syringe filters GF/PTFE-45/25 membrane material: glass fiber/polytetrafluoroethylene pore dia.: 1.0/0.45 μm , membrane dia.: 25 mm colorless PP case, labeled pack of 100



Code	Description	Packaging
MNL729271	CHROMAFIL Xtra GF/PTFE-45/25	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

