



Colonna Chromab. C18 ec, 3 mL, 500mg,ASP

Macherey Nagel

CHROMABOND columns C18 ec volume: 3 mL, content of sorbent: 500 mg
material: PP, with PE filter elements incl. special plug for ASPEC automate pack
of 100



Code	Description	Packaging
MNL730013ASP	Chromab. columns C18 ec, 3 mL, 500mg,ASP	100 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

