



Colonne Chromab. SiOH, 3 mL, 500 mg, ASP

Macherey Nagel

CHROMABOND columns SiOH volume: 3 mL, content of sorbent: 500 mg material: PP, with PE filter elements incl. special plug for ASPEC automat pack of 100





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

